

**d.) Remarks**

Any inquiry concerning this communication or earlier communications from the applicant should be directed to Chuan Li whose telephone number is (858) 361-7231. The applicant can normally be reached from 9:00 a.m. to 5:00 p.m. pacific standard time.

The applicant may also be reached at Expression Technologies Inc. at (858) 558-1861 or by fax at (858) 558-1883 or by email at [chuanli@exptec.com](mailto:chuanli@exptec.com).

Applicant Name: Chuan Li

Signature: 

Date: December 17, 2009

## BIO•SYNTHESIS

Lot No: B716-1

## Oligo Data Sheet

Date Created: 3/11/98  
Your Reference ID: OLIGO 1 1C2C01351  
Primer Lot Number: B716-1  
Author: MD  
Synthesis Scale: 50 nmole  
Primer Sequence (5' to 3'): CGC CCG CCG CCC GGG CGC CCC GCC TTC CGC  
TTC CTC GCT CAC TG

## Primer Data

Primer Length: 44  
Type: DNA  
Composition:      A              C              G              T              Others  
                         1              25              11              7              0  
                         2.3%          56.8%          25.0%          15.9%          0.0%

Molecular Weight (Ammonium Salt): 13231.8  
Exact Weight per OD (Ammonium Salt): 37.87  
Nanomoles per OD (Ammonium Salt): 2.86  
Micromolar Extinction Coefficient: 349.38  
Total ODs in This Tube: 5  
Total Amount in ug: 189.36  
Total Amount in nmoles: 14.31  
Purification: Desalted  
Melting Temperature in Celsius: 160.0

5' END OH  
3' END OH

Note: OD WILL VARY

biosyn® biosyn.com

800 DNA EXAM

Your source for custom DNA, peptides and molecular biology products

## BIO•SYNTHESIS

Lot No: B716-3

*Oligo Data Sheet*

Date Created: 3/11/98  
Your Reference ID: OLIGO 3 102001353  
Primer Lot Number: B716-3  
Author: MD  
Synthesis Scale: 50 nmole  
Primer Sequence (5' to 3'): CGC CCG CCG CCC GGG CGC CCC GCC AAC GCG  
GAA GTC AGC GCC CT

---

**Primer Data**

Primer Length:	44			
Type:	DNA			
Composition:				
A	C	G	T	Others
5	23	14	2	0
11.4%	52.3%	31.8%	4.5%	0.0%
Molecular Weight (Ammonium Salt):	13372.8			
Exact Weight per OD (Ammonium Salt):	34.85			
Nanomoles per OD (Ammonium Salt):	2.61			
Micromolar Extinction Coefficient:	383.67			
Total ODs in This Tube:	5			
Total Amount in ug:	174.27			
Total Amount in nmoles:	13.03			
Purification:	Desalted			
Melting Temperature in Celsius:	162.0			

---

5' END	OH
3' END	OH

Note: OD WILL VARY

*Your source for custom DNA, peptides and molecular biology products*

Sequence ID NO: 3

# BIO•SYNTHESIS

Lot No: B716-4

## Oligo Data Sheet

Date Created: 3/11/98  
Your Reference ID: OLIGO 4 102C 0331  
Primer Lot Number: B716-4  
Author: MD  
Synthesis Scale: 50 nmole  
Primer Sequence (5' to 3'): CGC CCG CCG CCC GGG CGC CCC GCC AAC GCA  
GAC CGT TCC GTG GC

### Primer Data

Primer Length: 44  
Type: DNA  
Composition:

A	C	G	T	Others
4	23	14	3	0
9.1%	52.3%	31.8%	6.8%	0.0%

Molecular Weight (Ammonium Salt): 13363.8  
Exact Weight per OD (Ammonium Salt): 35.38  
Nanomoles per OD (Ammonium Salt): 2.65  
Micromolar Extinction Coefficient: 377.73  
Total ODs in This Tube: 5  
Total Amount in ug: 176.9  
Total Amount in nmoles: 13.24  
Purification: Desalted  
Melting Temperature in Celsius: 162.0

5' END OH  
3' END OH

Note: OD WILL VARY

*Your source for custom DNA, peptides and molecular biology products*

# BIO•SYNTHESIS

## Oligo Data Sheet

Lot No: B716-5

**Date Created:** 3/11/98  
**Your Reference ID:** OLIGO 5 3C 015'  
**Primer Lot Number:** B716-5  
**Author:** MD  
**Synthesis Scale:** 50 nmole  
**Primer Sequence (5' to 3'):** CCG CCG CGC CGC TTC CAC TGA GCG TCA GAC CC

### Primer Data

**Primer Length:** 32  
**Type:** DNA  
**Composition:**

A	C	G	T	Others
4	16	8	4	0
12.5%	50.0%	25.0%	12.5%	0.0%

**Molecular Weight (Ammonium Salt):** 9668.4  
**Exact Weight per OD (Ammonium Salt):** 34.97  
**Nanomoles per OD (Ammonium Salt):** 3.62  
**Micromolar Extinction Coefficient:** 276.48  
**Total ODs in This Tube:** 5  
**Total Amount in ug:** 174.85  
**Total Amount in nmoles:** 18.08  
**Purification:** Desalted  
**Melting Temperature in Celsius:** 112.0

**5' END** OH  
**3' END** OH

**Note:** OD WILL VARY

*Your source for custom DNA, peptides and molecular biology products*

## BIO•SYNTHESIS

Lot No: B716-7

*Oligo Data Sheet*

Date Created: 3/11/98  
Your Reference ID: OLIGO 7 IGAN  
Primer Lot Number: B716-7  
Author: MD  
Synthesis Scale: 50 nmole  
Primer Sequence (5' to 3'): GGG CGG CGG GCG TTC GGG GAA ATG TGC GCG GA

---

**Primer Data**

Primer Length: 32  
Type: DNA  
Composition:

A	C	G	T	Others
4	6	18	4	0
12.5%	18.8%	56.3%	12.5%	0.0%

Molecular Weight (Ammonium Salt): 10068.4  
Exact Weight per OD (Ammonium Salt): 31.85  
Nanomoles per OD (Ammonium Salt): 3.16  
Micromolar Extinction Coefficient: 316.08  
Total ODs in This Tube: 5  
Total Amount in ug: 159.27  
Total Amount in nmoles: 15.82  
Purification: Desalted  
Melting Temperature in Celsius: 112.0

---

5' END OH  
3' END OH

Note: OD WILL VARY

*Your source for custom DNA, peptides and molecular biology products*

# BIO•SYNTHESIS

## Oligo Data Sheet

Lot No: B716-8

Date Created: 3/11/98  
Your Reference ID: OLIGO 8 1GKN  
Primer Lot Number: B716-8  
Author: MD  
Synthesis Scale: 50 nmole  
Primer Sequence (5' to 3'): GGG CGG CGG GCG TTG TCG GGA AGA TGC GTG  
AT

---

### Primer Data

Primer Length: 32  
Type: DNA  
Composition:

A	C	G	T	Others
4	5	17	6	0
12.5%	15.6%	53.1%	18.8%	0.0%

Molecular Weight (Ammonium Salt): 10058.4  
Exact Weight per OD (Ammonium Salt): 31.95  
Nanomoles per OD (Ammonium Salt): 3.18  
Micromolar Extinction Coefficient: 314.82  
Total ODs in This Tube: 5  
Total Amount in ug: 159.75  
Total Amount in nmoles: 15.88  
Purification: Desalted  
Melting Temperature in Celsius: 108.0

---

5' END OH  
3' END OH

Note: OD WILL VARY

*Your source for custom DNA, peptides and molecular biology products*

# BIO•SYNTHESIS

## Oligo Data Sheet

Lot No: B716-9

**Date Created:** 3/11/98  
**Your Reference ID:** OLIGO 9 14TN  
**Primer Lot Number:** B716-9  
**Author:** MD  
**Synthesis Scale:** 50 nmole  
**Primer Sequence (5' to 3'):** GGG CGG CGG GCG TTC TCA TGT TTG ACA GCT TA

### Primer Data

**Primer Length:** 32  
**Type:** DNA  
**Composition:**

A	C	G	T	Others
4	7	12	9	0
12.5%	21.9%	37.5%	28.1%	0.0%

**Molecular Weight (Ammonium Salt):** 9903.4  
**Exact Weight per OD (Ammonium Salt):** 33.11  
**Nanomoles per OD (Ammonium Salt):** 3.34  
**Micromolar Extinction Coefficient:** 299.07  
**Total ODs in This Tube:** 5  
**Total Amount in ug:** 165.57  
**Total Amount in nmoles:** 16.72  
**Purification:** Desalted  
**Melting Temperature in Celsius:** 102.0

**5' END** OH  
**3' END** OH

**Note:** OD WILL VARY

*Your source for custom DNA, peptides and molecular biology products*



## BIO•SYNTHESIS

Lot No: B716-10

*Oligo Data Sheet*

Date Created: 3/11/98  
Your Reference ID: OLIGO 10  
Primer Lot Number: B716-10  
Author: MD  
Synthesis Scale: 50 nmole  
Primer Sequence (5' to 3'): GGG CGG CGG GCG AAG CCA CTG GAG CAC CTC AA

---

**Primer Data**

Primer Length:	32				
Type:	DNA				
Composition:	A	C	G	T	Others
	7	10	13	2	0
	21.9%	31.3%	40.6%	6.3%	0.0%
Molecular Weight (Ammonium Salt):	9910.4				
Exact Weight per OD (Ammonium Salt):	31.42				
Nanomoles per OD (Ammonium Salt):	3.17				
Micromolar Extinction Coefficient:	315.45				
Total ODs in This Tube:	5				
Total Amount in ug:	157.08				
Total Amount in nmoles:	15.85				
Purification:	Desalted				
Melting Temperature in Celsius:	110.0				

---

5' END	OH
3' END	OH

Note: OD WILL VARY

*Your source for custom DNA, peptides and molecular biology products*

# BIO•SYNTHESIS

## Oligo Data Sheet

Lot No: B716-11

**Date Created:** 3/11/98  
**Your Reference ID:** OLIGO 11 3GAC  
**Primer Lot Number:** B716-11  
**Author:** MD  
**Synthesis Scale:** 50 nmole  
**Primer Sequence (5' to 3'):** GCG GCG CGG CGG TAC GGG GTC TGA CGC TCA GT

### Primer Data

**Primer Length:** 32  
**Type:** DNA  
**Composition:**

A	C	G	T	Others
3	9	15	5	0
9.4%	28.1%	46.9%	15.6%	0.0%

**Molecular Weight (Ammonium Salt):** 9939.4  
**Exact Weight per OD (Ammonium Salt):** 33.32  
**Nanomoles per OD (Ammonium Salt):** 3.35  
**Micromolar Extinction Coefficient:** 298.26  
**Total ODs in This Tube:** 5  
**Total Amount in ug:** 166.62  
**Total Amount in nmoles:** 16.76  
**Purification:** Desalted  
**Melting Temperature in Celsius:** 112.0

**5' END** OH  
**3' END** OH

**Note:** OD WILL VARY

*Your source for custom DNA, peptides and molecular biology products*

Sequence ID NO: 10

# BIO•SYNTHESIS

## Oligo Data Sheet

Lot No: B716-12

Date Created: 3/11/98  
Your Reference ID: OLIGO 12 3GKC  
Primer Lot Number: B716-12  
Author: MD  
Synthesis Scale: 50 nmole  
Primer Sequence (5' to 3'): GCG GCG CGG CGG ATC GCC CCA TCA TCC AGC CA

### Primer Data

Primer Length: 32  
Type: DNA  
Composition:

A	C	G	T	Others
5	14	10	3	0
15.6%	43.8%	31.3%	9.4%	0.0%

Molecular Weight (Ammonium Salt): 9757.4  
Exact Weight per OD (Ammonium Salt): 33.61  
Nanomoles per OD (Ammonium Salt): 3.44  
Micromolar Extinction Coefficient: 290.34  
Total ODs in This Tube: 5  
Total Amount in ug: 168.03  
Total Amount in nmoles: 17.22  
Purification: Desalted  
Melting Temperature in Celsius: 112.0

5' END OH  
3' END OH

Note: OD WILL VARY

*Your source for custom DNA, peptides and molecular biology products*

**BIO•SYNTHESIS**  
*Oligo Data Sheet*

Lot No: B716-13

**Date Created:** 3/11/98  
**Your Reference ID:** OLIGO 13 39TC  
**Primer Lot Number:** B716-13  
**Author:** MD  
**Synthesis Scale:** 50 nmole  
**Primer Sequence (5' to 3'):** GCG GCG CGG CGG TTC ACG TTC GCT CGC GTA TC

---

**Primer Data**

**Primer Length:** 32  
**Type:** DNA  
**Composition:**

A	C	G	T	Others
2	11	12	7	0
6.3%	34.4%	37.5%	21.9%	0.0%

**Molecular Weight (Ammonium Salt):** 9825.4  
**Exact Weight per OD (Ammonium Salt):** 34.87  
**Nanomoles per OD (Ammonium Salt):** 3.55  
**Micromolar Extinction Coefficient:** 281.79  
**Total ODs in This Tube:** 5  
**Total Amount in ug:** 174.34  
**Total Amount in nmoles:** 17.74  
**Purification:** Desalted  
**Melting Temperature in Celsius:** 110.0

---

**5' END** OH  
**3' END** OH

**Note:** OD WILL VARY

*Your source for custom DNA, peptides and molecular biology products*

## BIO•SYNTHESIS

Lot No: B716-14

*Oligo Data Sheet*

**Date Created:** 3/11/98  
**Your Reference ID:** OLIGO 14 3GCC  
**Primer Lot Number:** B716-14  
**Author:** AH  
**Synthesis Scale:** 50 nmol  
**Primer Sequence (5' to 3'):** GCG GCG CGG CGG AAG CAC ACG GTC ACA CTG CT

**Primer Data**

**Primer Length:** 32  
**Type:** DNA  
**Composition:**

A	C	G	T	Others
6	11	12	3	0
18.8%	34.4%	37.5%	9.4%	0.0%

**Molecular Weight (Ammonium Salt):** 9861.4  
**Exact Weight per OD (Ammonium Salt):** 32.27  
**Nanomoles per OD (Ammonium Salt):** 3.27  
**Millimolar Extinction Coefficient:** 305.55  
**Total ODs in This Tube:** 5  
**Total Amount in ug:** 161.37  
**Total Amount in nmoles:** 16.36  
**Purification:** Desalted  
**Melting Temperature in Celsius:** 110.0

**5' END** OH  
**3' END** OH

**Note:** OD WILL VARY

*Your source for custom DNA, peptides and molecular biology products*

↓ 14K 10<sup>5</sup>

↓ 1ml 70% EtOH

↓ 2<sup>1</sup>↓ air dry  
pellets

The Amount of DNA of #2 and #4 products is about 18 µg.  
 amount of DNA of #7, #8, #9 and #4 products is about 28 µg.  
 90 µl and 140 µl H<sub>2</sub>O to make their final concentration at about 0.2 µg/µl.  
 p them at -20°C.

4 PCR OS

6-19-98

Now I am doing PCR analysis since the PCR Test worked fine.

DNA	Primers	Product	bps	Comments
pBR322, PstI	3+5	01	1535	
"	1+5	02	813	
pUC19, RI	1+5	03	813	
pACYC177, BHI	4+5	04	1011	
"	1+5	05	726	
pBR322, PstI	2+5	06	1057	
pACYC177, BHI	1+6	07	694	
"	4+6	08	979	
pACYC184, BHI	1+6	09	694	
pBR322, PvuII	7+11	S1	1130	
pUC19, RI	7+11	S2	1130	
pACYC177, BHI	7+11	S3	1130	
"	8+12	S4	1219	
pBR322, PvuII	9+13	S5	1552	
pACYC184, BHI	10+14	S6	1104	

DNA are from page 6. They are diluted to 1 ng/µl and 1 µl was  
 1 for each reaction.

Primers are from page OS. The estimated concentration is about  
 8 µl. 2 µl of each primer was used.

Dilute the DNA into appropriate concentration.

Making master mix as on page 6 except the number of reactions  
 1d be 16 instead of 12. One extra is for negative To Page No. 9

sed & Understood by me,

Date

6/18/98

Invented by

Recorded by


Date

6-19-98

Page No. 8 control.

The PCR condition is same as on page 7.

100ul of each PCR product



6/24/98  
↓ take out 3ul

↓ run on 0.8% agarose TBE


← ↓ take a picture

All the reactions worked. However, some worked better than others.

#1 to 5 and 10 to 14 seems have many non-specific products. By increasing the template concentration, specific product may be increased and non-specific products may be decreased.

6-20-98 → Repeat reactions #1 to 5 and 10 to 14 by increasing the template 10X (use 10ng)

100ul of each PCR product.



6/24/98  
↓ take out 3ul

↓ run on 0.8% agarose TBE

← ↓ take a picture

All the reactions appear to be better except #13. I probably made mistake when adding template or primers for #13. Keep the PCR products at -20°C

## 5. PCR SI

6-21-98

DNA	Amount	Primers	P	bps	Comments
pACYC117 BHI	1ng	8+12	S4	1219	
"	10ng	"	"	"	
pACYC117 pAET	1ng	"	"	"	
"	10ng	"	"	"	
pGEX-3X	1ng	15+17	I	1198	

The primers are from page 5. 2ul of each primer was used. The PCR condition is same as on page 7. To Page No. 10

Used & Understood by me,

*Steven F. Gessert*  
STEVEN F. GESSERT

Date

6/24/98

Invented by

*Chuan Li*  
Recorded by *Chuan Li*

Date

6-20-98  
6-21-98

Page No. 92

(2) By comparing the yields of the minipreps ~~on~~ between page 11 and 92, it is clear that the copy number of plasmids is determined by the selection marker it contains. In OS4 minipreps, the A. based plasmids give very low yields except #15.

(3) ~~The~~ SmaI has two sites on OS4 constructs. Therefore two bands generated after SmaI digestion. However the size of the DNA appear to be different even from the colonies ~~are~~ picked from same plate. This observation may be artifacts of electrophoresis, but I need pay attention on this observation in further analysis.

## 22. OS Mediprep Test 1 9.8.98

Use the residual O/N cultures from 4 ml inoculation (these residual ones are kept at 4°C) to seed 5 ml LB with appropriate antibiotic - prep number, page and O/N culture number are indicated in the following table.

1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
84	84	<del>84</del> 87	87	85	87	85		91	<del>88</del> 81	91			92		
2	28	1	24	18	22	21	25	19	8	13	16	2	8	10	15
Amp				Tet				chl				Kan			

The above table is messy. I will re-prepare the table below:

1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
84	84	<del>84</del> 87	87	85	87	85		91		87	91			92	
2	28	16	4	18	22	21	25	19	8	13	16	2	8	10	15
Amp				Tet				chl				Kan			

5 ml inoculated LB w/ appropriate antibiotic

↓ 37°C w/ shaking at 300 RPM for 3h,

1 to 4, 5, 7, 13 and 15 grow fast.

6, 8, 9 to 12, 14 and 16 grow slowly.

↓ 37°C for another 3h w/ shaking

9 to 12, 14 and 16 still do not grow well

↓ 37°C w/ shaking O/N. ~13hs

9.9.98

inoculate all 5 ml into 50 ml LB w/ antibiotic

↓ 37°C for 1h w/ shaking

Test OD590

#1 : 0.435

#12 : 0.278

To Page No. 94

Used &amp; Understood by me,

Date

Invented by

Date

Recorded by

9.8.98



Page No. 93

↓ 37°C for another 1 h w/ shaking.

TEST OD<sub>575</sub>

#1: 0.718

#12: 0.283

#8: 0.876

#16: 0.257

the final concentration of chl  
 $5 \times 35 / 50 \approx 10 \mu\text{g/ml}$

↓ add 15  $\mu\text{l}$  35 mg/ml chl to 1-8 and 13-16

↓ 37°C w/ shaking for 4 hrs

#1: 1.260

#12: 1.192

#8: 0.992

#16: 0.337

↓ 3.2 K for 15' (actually 10' should be enough)

pellets (#16 has smallest pellet)

↓ 0.5 ml LETT

↓ completely resuspend cells by pipetting

↓ transfer to eppendorf tubes

↓ boil for 90" → 120"

↓ 14K 10'

\*1 → \*2 → \*4

↓ transfer supernatants to new tubes

supernatants \*3

↓ 1V  $\phi$ -chl ext.

↓ 2V EtOH

↓ 5' 14K

pellets

↓ 70% EtOH wash

pellets

↓

also has a small pellet.  
 and 11 have biggest pellets,  
 are most viscous after  
 resuspend the cells

after 10' spinning, #1-5, 7, 8,  
 13-15 did not pellet, boil  
 for another 3' and repeat  
 spinning.

after second spin, they form  
 a pellet. Freeze them in  
 ice/EtOH bath, thaw them,  
 spin for another 10'  
 #5, 7, 14, 15 are treated  
 dry ice/EtOH bath)

adjust the supernatants to  
 same volume by LETT

To Page No. 95

Read &amp; Understood by me,

Date

Invented by

Date

Recorded by

9.9.98

om Page No. 94 \* 4 The lysed cell pellets are different in size. The approximate pellet size (CAPS) are listed in the following table:

cap#	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15
PS	400	900	300	250	400	100	250	200	200	200	150	120	100	250	250

The approximate pellet sizes (CAPS) are in microliters.

pellets after 70% EtOH wash are different in size. Their relative size are listed in the following table:

q#	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15
re	L	L	L	L	L	<del>M</del> S	L	M	S	S	S	VS	L	L	L

↓ air dry the pellets O/N.

↓ resuspend in 200ul TE w/RNase A 9.10

After O/N air dry, the pellets are difficult to be dissolved, especially the larger pellets. When they are finally dissolved (take about 2 hrs w/ vortexing), they form heavy foams while vortexing.

Add another 200ul TE w/RNase A to large pellet tubes namely # 1-5, 7, 13-15.

It is amazing that all the pellets seem dissolve completely these efforts.

Sma I digestion:

Master Mix for each RXN	16 RXNs
H <sub>2</sub> O 6.9 ul	110.4
10X Buffer 1 ul	16
Sma I 0.1 ul	1.6
8 ul/RXN	128 ul totally


take out 2ul from each sample

↓ add 8ul Master mix

↓ R.T for 1h

↓ run on 0.8% agarose TBE

↓ take a picture. (416, 1/2 sec,



Result analysis:

- ① Genomic DNA contamination is serious possible solutions
- ② Decrease the boiling time
- ③ Do not use pipette to resuspend the cell
- ④ Decrease the Triton-X100 Concentration.
- ⑤ Most of the preps have enough DNA for future usage. However 3, 7 and 15 appear to have very little DNA (0.025 per band).

possible solutions ⑥ Use higher concentration of chloroform/phenol on To Page No.

Informed & Understood by me,

Date

Invented by

Date

Recorded by

9-9-10-98