

Page No. 7

↓ 14K 10  
↓ 1ml 70% EtOH  
↓ 2'  
↓ air dry  
pellets

The Amount of DNA of #2 and #4 products is about 18 ug.  
Amount of DNA of #7, #8, #9 and #4 products is about 28 ug.  
90ul and 140ul H<sub>2</sub>O to make their final concentration at about 0.2 ug/ul.  
p them at -20°C.

4. PCR OS

6-19-98

Now I am doing PCR analysis since the PCR Test worked fine.

DNA	Primers	Product	bps	Comments
pBR322, PstI	3+5	01	1535	
"	1+5	02	813	
pUC19, RI	1+5	03	813	
pACYC177, BHI	4+5	04	1011	
"	1+5	05	726	
pBR322, PstI	2+5	06	1057	
pACYC177, BHI	1+6	07	694	
"	4+6	08	979	
pACYC184, BHI	1+6	09	694	
pBR322, PvuII	7+11	S1	1130	
pUC19, RI	7+11	S2	1130	
pACYC177, BHI	7+11	S3	1130	
"	8+12	S4	1219	
pBR322, PvuII	9+13	S5	1552	
pACYC184, BHI	10+14	S6	1104	

DNA are from page 6. They are diluted to 1 ng/ul and 1 ul was  
1 for each reaction.

Primers are from page 15. The estimated concentration is about  
g/ul. 2 ul of each primer was used.

Dilute the DNA into appropriate concentration.

Making master mix as on page 6. except the number of reactions  
ld be 16 instead of 12. One extra is for negative. To Page No. 9

sed & Understood by me,

*[Signature]*

Date

6/18/98

Invented by

*[Signature]*

Date


6-19-98

Recorded by

*[Signature]*

Page No. 8 control.

The PCR condition is same as on page 7.


100  $\mu$ l of each PCR product6/24/98 → take out 3  $\mu$ l

6/24/98 → run on 0.8% agarose TBE

← take a picture

All the reactions worked. However, some worked better than others.

#1 to 5 and 10 to 14 seems have many non-specific products. By increasing the template concentration, specific product may be increased and non-specific products may be decreased.

6-20-98 → Repeat reactions #1 to 5 and 10 to 14 by raising the template 10x (use 10  $\mu$ l)100  $\mu$ l of each PCR product.6/24/98 → take out 3  $\mu$ l

6/24/98 → run on 0.8% agarose TBE

← take a picture

All the reactions appear to be better except #13. I probably made mistake when adding template or primers for #13. Keep the PCR products at -20°C

## 5. PCR SI

6-21-98

DNA	Amount	Primers	P	bps	Comments
PACYC177 BH1	1 $\mu$ g	8+12	S4	1219	
"	10 $\mu$ g	"	"	"	
PACYC177 PAEII	1 $\mu$ g	"	"	"	
"	10 $\mu$ g	"	"	"	
PGEX -3X	1 $\mu$ g	15+17	I	1198	

The primers are from page 5. 2  $\mu$ l of each primer was used.

The PCR condition is same as on page 7.

To Page No. 10

Used &amp; Understood by me,

STEVEN F. GESSERT

Date

6/24/98

Invented by

CHUAN LI

Recorded by

CHUAN LI

Date

6-20-98

6-21-98

Page No. 92

② By comparing the yields of the minipreps between page 11 and 92, it is clear that the copy number of plasmids is determined by the selection marker it contains. In OS4 minipreps, the A. based plasmids give very low yields except #15.

③ ~~The~~ SmaI has two sites on OS4 constructs. Therefore two bands generated after SmaI digestion. However the size of the DNA appear be different even from the colonies ~~are~~ picked from same plate. This observation may be artifacts of electrophoresis, but I need pay attention on this observation in further analysis.

## 22. OS Mediprep Test 1 9.8.98

Use the residual O/N cultures from 4 ml inoculation (these residual ones are kept at 4°C) to seed 5 ml LB with appropriate antibiotic prep number, page and O/N culture number are indicated in the following table.

1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
84	84	<del>84</del> 87	87	85	87	85		91	<del>84</del> 88	<del>84</del> 83	91		92		
2	28	1	24	18	22	21	25	19	8	13	16	2	8	10	15
Amp				Tet				chl				Kan			

The above table is messy. I will re-prepare the table below:

1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
84	84	<del>84</del> 87	87	85	87	85		91		87	91			92	
2	28	16	4	18	22	21	25	19	8	13	16	2	8	10	15
Amp				Tet				chl				Kan			

5 ml inoculated LB w/ appropriate antibiotic

↓ 37°C w/ shaking at 300 RPM for 3h,

1 to 4, 5, 7, 13 and 15 grows fast.

6, 8, 9 to 12, 14 and 16 grow slowly.

↓ 37°C for another 3h w/ shaking

9 to 12, 14 and 16 still do not grow well

↓ 37°C w/ shaking O/N. ~13hs

9.9.98

inoculate all 5 ml into 50 ml LB w/ antibiotic

↓ 37°C for 1h w/ shaking

Test OD590

#1 : 0.435

#12 : 0.278

To Page No. 94

Used & Understood by me,

Date

Invented by

Date

Recorded by

9-11-98

Chuan Li

9.8.98

Page No. 93

↓ 37°C for another 1h w/ shaking.

Test OD<sub>525</sub>

#1: 0.718      #12: 0.283  
 #8: 0.876      #16: 0.257

the final concentration of chl  
 5 x 35/50 ≈ 10 ug/ml

↓ add 15 ul 35 mg/ml chl to 1-8 and 13-16

↓ 37°C w/ shaking for 4hs

#1: 1.260      #12: 1.192  
 #8: 0.992      #16: 0.337

↓ 32K for 5' (actually 10' should be enough)

also has a small pellet.  
 and 11 have biggest pellets,  
 are most viscous after  
 suspend the cells

pellets (#16 has smallest pellet)

↓ 0.5 ml LETT

↓ completely resuspend cells by pipetting

↓ transfer to eppendorf tubes

↓ boil for 10' → 120"

after 10' spinning, #1-5, 7, 8,  
 13-15 did not pellet, boil  
 for another 3' and repeat  
 spinning.

↓ 14K 10'

\*1 → \*2 → \*4

↓ transfer supernatants to new tubes

supernatants \*3

after second spin, they form  
 a pellet. Freeze them in  
 ice/EtOH bath, thaw them,  
 spin for another 10'

↓ 1V φ = chl ext.

↓ 2V EtOH

\*5, 7, 14, 15 are treated  
 dry ice/EtOH bath)

↓ 5' 14K

adjust the supernatants to  
 same volume by LETT

pellets

↓ 70% EtOH wash

pellets

↓

To Page No. 95

Read &amp; Understood by me,

Date

Invented by

Date

Recorded by

9-11-98

Chuan Li

9-9-98

Chuan Li

om Page No. 94 \*4 The lysed cell pellets are different in size. The approximate pellet size (CAPS) are listed in the following table:

cap #	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15
PS	400	900	300	250	400	100	250	200	200	200	150	120	100	250	250

The approximate pellet sizes (CAPS) are in microliters.

pellets after 70% EtOH wash are different in size. Their relative size are listed in the following table:

#	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15
rel	L	L	L	L	L	<del>MS</del> S	L	M	S	S	S	VS	L	L	L

↓ air dry the pellets o/n.

↓ resuspend in 200ul TE w/RNase A 9:10.

After o/n air dry, the pellets are difficult to be dissolved especially the larger pellets. When they are finally dissolved (take about 2 hrs w/ vortexing), they form heavy foams while vortexing.

Add another 200ul TE w/RNase A to large pellet tubes namely #1-5, 7, 13-15.

It is amazing that all the pellets seem dissolve completely w/ these efforts.

Sma I digestion:

Master Mix for each RXN	16 RXNs
H <sub>2</sub> O 6.9 ul	110.4
10X Buffer 1 ul	16
Sma I 0.1 ul	1.6
8 ul/RXN	128 ul totally.


take out 2ul from each sample

↓ add 8ul master mix

↓ R.T for 1h

↓ run on 0.8% agarose TBE

↓ take a picture. (716, 1/2 sec)



Result analysis:

- ① Genomic DNA contamination is serious possible solutions
- ② Decrease the boiling time
- ③ Do not use pipette to resuspend the cell
- ④ Decrease the Triton-X100 Concentration.
- ⑤ Most of the preps have enough DNA for future usage. (However 3, 7 and 15 appear to have very little DNA GO-0.25 per band).

possible solutions ⑥ Use higher concentration of chloroform on

Witnessed & Understood by me,

Date

Invented by

Date

Recorded by

9-9-10-98

To Page No. 9